

TRIOS Core Wired

Core. Essential functionality.

- ✓ High accuracy and precision
- ✓ Realistic scan colors
- ✓ Automatic scan clean-up
- ✓ Print-ready model design
- ✓ Unrivalled open dental-ecosystem
- ✓ Opportunity for in-house design and production

TECHNICAL SPECIFICATIONS

GENERAL

What's in the box*	TRIOS Core scanner wand, 5x TRIOS 3 scanner tips, TRIOS 3 protection tip, TRIOS 3 calibration tip, TRIOS 3 color calibration tip, TRIOS 3 color calibration kit, TRIOS 3 scanner pod, 3Shape USB-C custom cable, 3Shape license dongle, USB dongle with software, documentation
Packaging Dimensions (LxWxH cm)	44 x 12,5 x 31
Invisalign Connection	No

HARDWARE

Scanner Weight	340 g
Scanner Dimensions (LxWxH cm):	27,3 x 4 x 4,8
Tip Dimensions (LxWxH cm):	8,6 x 2 x 2,05
Scanner Cord Length (cm)	150
Pod Dimensions (DxH cm)	12 x 7
Pod Weight	510 g
USB-C Cable Length (cm)	100

* TRIOS Core is compatible with TRIOS 3 accessories



3shape 

Let's change dentistry together