

myQuickmat Forte starter kit

new!

QuickmatFLEX
titanium sectional
matrices included

Advantages

Restorative essentials for
the recreation of tight,
anatomically-correct contact
points, and natural contours

Selected restorative solutions
for improved matrix retention
and adaptation

Great choice of sectional
matrices to optimise your
Class II restorations

One system for restoring teeth
with different crown height
and wide cavity preps

sectional
matrix
systems



polydentia
swiss manufacture 



myQuickmat Forte starter kit

power up your
restorations

Testimonial



myQuickmat Forte sectional matrix system provides top-class solutions for Class II restorations. The separator ring myRing Forte provides effective teeth separation and matrix retention. In combination with myWedge, the system ensures a perfect matrix fit and reproduction of tight, anatomical contact points. The addition of the ultra-thin and resilient QuickmatFLEX sectional matrices further enhances the versatility of this system.


Dr. Giuseppe Chiodera,
Italy



Images courtesy by Dr Osama Shaalan



for more
information



new! QuickmatFLEX

Titanium sectional matrices with advanced elastic memory. Ultra-thin for easy insertion into tight proximal spaces and resilient for optimal deformation resistance



Highly resistant
to deformation

Ultra-thin
titanium
(0.03 mm)

myWedge

V-shaped wedge made of hi-tech plastic material. Perfectly compresses and adapts to the interproximal anatomy of teeth



LumiContrast

Dark-blue steel matrices which significantly increase contrast and reduce glare effect, especially while using loupes or microscopes

myTines

Replaceable, interchangeable and autoclavable plastic extremities. myTines make myRing Forte adaptable to teeth with different crown heights as well as teeth with large preparations or missing cusps.



Obtain anatomically correct Class II cavity restorations with tight contact points and natural contours using the restorative solutions included in myQuickmat Forte starter kit.

The system allows you to restore teeth with different crown heights as well as teeth with large cavity preparations or missing cusps thanks to the replaceable extremity system myTines.

All components of the kit have been selected to optimise the restorative procedure, allowing quick and predictable results with minimal finishing steps.

myQuickmat Forte starter kit - REF 6883

- 2 pcs. myRing Forte, 2 pairs myTines Small (S),
- 1 pair myTines Medium (M),
- 1 pair myTines Large (L),
- 10 pcs. QuickmatFLEX Sectional Matrices, Premolar 0.03mm / 5mm,
- 10 pcs. QuickmatFLEX Sectional Matrices, Molar 0.03mm / 6.4mm,
- 5 pcs. QuickmatFLEX Sectional Matrices, Molar Deep 0.03mm / 6.4mm + Ext.,
- 10 pcs. LumiContrast Sectional Matrices, Premolar 0.04mm / 5mm,
- 10 pcs. LumiContrast Sectional Matrices, Molar 0.04mm / 6.4mm,
- 5 pcs. LumiContrast Sectional Matrices, Molar Deep 0.04mm / 6.4mm + Ext.,
- 10 pcs. myWedge, Extrasmall (XS),
- 10 pcs. myWedge, Small (S),
- 10 pcs. myWedge, Medium (M),
- 10 pcs. myWedge, Large (L)

